

SEC16B Antibody



PACO39850

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Required for secretory cargo traffic from the endoplasmic reticulum to the Golgi apparatus and for normal transitional endoplasmic reticulum (tER) organization.

Gene ID:

SEC16B

Uniprot

Q96JE7

Synonyms:

Protein transport protein Sec16B (Leucine zipper transcription regulator 2) (Regucalcin gene promoter region-related protein p117) (RGPR-p117) (SEC16 homolog B), SEC16B, KIAA1928 LZTR2 RGPR SEC16S


Immunogen:

Recombinant Human Protein transport protein Sec16B protein (5-304AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39850 at dilution of 1:100.