

FAM71F1 Antibody




PACO39854

Product Information

Size:	Protein Background:
50ug	nucleus.
Reactivity:	Gene ID:
Human	FAM71F1
Source:	Uniprot
Rabbit	Q96KD3
Isotype:	Synonyms:
IgG	Protein FAM71F1 (Protein FAM137A) (Testis development protein NYD-SP18), FAM71F1, FAM137A
Applications:	Immunogen:
ELISA, IHC	Recombinant Human Protein FAM71F1 protein (121-343AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO39854 at dilution of 1:100.