

PNMA5 Antibody



PACO39882

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

PNMA5

Uniprot

Q96PV4

Synonyms:

Paraneoplastic antigen-like protein 5 (Tumor antigen BJ-HCC-25), PNMA5, KIAA1934


Immunogen:

Recombinant Human Paraneoplastic antigen-like protein 5 protein (4-442AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO39882 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO39882 at dilution of 1:100.