

LILRA2 Antibody, FITC conjugated



PACO39896

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Does not bind class I MHC antigens.

Gene ID:

LILRA2

Uniprot

Q8N149

Synonyms:

Leukocyte immunoglobulin-like receptor subfamily A member 2 (CD85 antigen-like family member H) (Immunoglobulin-like transcript 1) (ILT-1) (Leukocyte immunoglobulin-like receptor 7) (LIR-7) (CD antigen CD85h), LILRA2, ILT1 LIR7

Immunogen:

Recombinant Human Leukocyte immunoglobulin-like receptor subfamily A member 2 protein (9-275AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A