

FLAD1 Antibody



PACO39954

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000,
IHC:1:20-1:200

Protein Background:

Catalyzes the adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme.

Gene ID:

FLAD1

Uniprot

Q8NFF5

Synonyms:

FAD synthase (EC 2.7.7.2) (FAD pyrophosphorylase) (FMN adenylyltransferase) (Flavin adenine dinucleotide synthase) [Includes: Molybdenum cofactor biosynthesis protein-like region; FAD synthase region], FLAD1


Immunogen:

Recombinant Human FAD synthase protein (1-490AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: FLAD1 antibody at 4 μ g/ml. Lane 1: K562 whole cell lysate. Lane 2: HepG2 whole cell lysate. Lane 3: MDA-MB-231 whole cell lysate. Lane 4: Hela whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 66, 55, 50, 33, 37 kDa. Observed band size: 66 kDa.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO39954 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO39954 at dilution of 1:100.