

DUSP19 Antibody



PACO39990

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Has a dual specificity toward Ser/Thr and Tyr-containing proteins.

Gene ID:

DUSP19

Uniprot

Q8WTR2

Synonyms:

Dual specificity protein phosphatase 19 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase TS-DSP1) (Low molecular weight dual specificity phosphatase 3) (LMW-DSP3) (Protein phosphatase SKRP1) (Stress-activated protein kinase pathway-regulating phosphatase 1) (SAPK pathway-regulating phosphatase 1), DUSP19, DUSP17, LMWDSP3, SKRP1


Immunogen:

Recombinant Human Dual specificity protein phosphatase 19 protein (1-217AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO39990 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO39990 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO39990 at dilution of 1:100.