

MAIP1 Antibody, Biotin conjugated



PACO40021

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MAIP1

Uniprot

Q8WWC4

Synonyms:

m-AAA protease-interacting protein 1, mitochondrial (Matrix AAA peptidase-interacting protein 1), MAIP1, C2orf47

Immunogen:

Recombinant Human m-AAA protease-interacting protein 1, mitochondrial protein (97-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A