

FAM3C Antibody



PACO40046

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May be involved in retinal laminar formation. Promotes epithelial to mesenchymal transition.

Gene ID:

FAM3C

Uniprot

Q92520

Synonyms:

Protein FAM3C (Interleukin-like EMT inducer), FAM3C, ILEI


Immunogen:

Recombinant Human Protein FAM3C protein (31-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: FAM3C antibody at 4 μ g/ml + A431 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 25 kDa. Observed band size: 25 kDa.