

SCAMP4 Antibody



PACO40066

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Probably involved in membrane protein trafficking.

Gene ID:

SCAMP4

Uniprot

Q969E2

Synonyms:

Secretory carrier-associated membrane protein 4 (Secretory carrier membrane protein 4), SCAMP4


Immunogen:

Recombinant Human Secretory carrier-associated membrane protein 4 protein (170-228AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO40066 at dilution of 1:100.



Immunofluorescent analysis of U251 cells using PACO40066 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).