

ROPN1L Antibody



PACO40094

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, motile cilium, epithelial cilium movement.

Gene ID:

ROPN1L

Uniprot

Q96C74

Synonyms:

Ropporin-1-like protein (ROPN1-like protein) (AKAP-associated sperm protein), ROPN1L, ASP


Immunogen:

Recombinant Human Ropporin-1-like protein (1-230AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO40094 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human breast cancer using PACO40094 at dilution of 1:100.