

# TRIM11 Antibody



PACO40126


---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	TRIM11
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q96F44
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	E3 ubiquitin-protein ligase TRIM11 (EC 2.3.2.27) (Protein BIA1) (RING finger protein 92) (RING-type E3 ubiquitin transferase TRIM11) (Tripartite motif-containing protein 11), TRIM11, RNF92
<b>Applications:</b>	<b>Immunogen:</b>
ELISA, IHC	Recombinant Human E3 ubiquitin-protein ligase TRIM11 protein (267-468AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
ELISA:1:2000-1:10000, IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human heart tissue using PACO40126 at dilution of 1:100.