

# NARS2 Antibody



PACO40138

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

mitochondrion, asparagine-tRNA ligase activity, asparaginyl-tRNA aminoacylation, gene expression, tRNA aminoacylation for protein translation.

**Gene ID:**

NARS2

**Uniprot**

Q96I59

**Synonyms:**

Probable asparagine--tRNA ligase, mitochondrial (EC 6.1.1.22) (Asparaginyl-tRNA synthetase) (AsnRS), NARS2

**Immunogen:**


Recombinant Human Probable asparagine--tRNA ligase, mitochondrial protein (15-210AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO40138 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO40138 at dilution of 1:100.