

METTL2A Antibody, Biotin conjugated



PACO40149

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable S-adenosyl-L-methionine-dependent methyltransferase that mediates 3-methylcytidine modification of some tRNAs.

Gene ID:

METTL2A

Uniprot

Q96IZ6

Synonyms:

Methyltransferase-like protein 2A (EC 2.1.1. -), METTL2A, METTL2

Immunogen:

Recombinant Human Methyltransferase-like protein 2A protein (1-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A