

ZHX1-C8orf76 Antibody



PACO40158

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

ZHX1-C8orf76

Uniprot

Q96EF9

Synonyms:

Zinc fingers and homeoboxes protein 1, isoform 2 (ZHX1-C8orf76 readthrough transcript protein), ZHX1-C8orf76, C8orf76


Immunogen:

Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 protein (1-292AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: ZHX1-C8orf76 antibody at 2µg/ml. Lane 1: Mouse heart tissue. Lane 2: Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 33 kDa. Observed band size: 33 kDa.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO40158 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO40158 at dilution of 1:100.