

ZHX1-C8orf76 Antibody, FITC conjugated



PACO40160

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

ZHX1-C8orf76

Uniprot

Q96EF9

Synonyms:

Zinc fingers and homeoboxes protein 1, isoform 2 (ZHX1-C8orf76 readthrough transcript protein), ZHX1-C8orf76, C8orf76

Immunogen:

Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 protein (1-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A