
Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Hydrolyzes a variety of dipeptides including L-carnosine but has a strong preference for Cys-Gly. Isoform 2 may play a role as tumor suppressor in hepatocellular carcinoma (HCC) cells. Ref.9 Ref.10.

Gene ID:

CNDP2

Uniprot

Q96KP4

Synonyms:

Cytosolic non-specific dipeptidase (EC 3.4.13.18) (CNDP dipeptidase 2) (Carnosine dipeptidase II) (Epididymis secretory protein Li 13) (Glutamate carboxypeptidase-like protein 1) (Peptidase A), CNDP2, CN2 CPGL HEL-S-13 PEPA


Immunogen:

Recombinant Human Cytosolic non-specific dipeptidase protein (2-475AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO40170 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO40170 at dilution of 1:100.