

PRMT6 Antibody



PACO40174

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:**Gene ID:**

PRMT6

Uniprot

Q96LA8

Synonyms:

Protein arginine N-methyltransferase 6 (EC 2.1.1.319) (Heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6) (Histone-arginine N-methyltransferase PRMT6), PRMT6, HRMT1L6


Immunogen:

Recombinant Human Protein arginine N-methyltransferase 6 protein (119-375AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: PRMT6 antibody at 0.4 μ g/ml. Lane 1: 293T whole cell lysate. Lane 2: MCF-7 whole cell lysate. Lane 3: A431 whole cell lysate. Lane 4: K562 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42, 33 kDa. Observed band size: 42 kDa.



Immunofluorescent analysis of HeLa cells using PACO40174 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO40174 at dilution of 1:100.