

CXorf58 Antibody



PACO40178

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CXorf58

Uniprot

Q96LI9

Synonyms:

Putative uncharacterized protein CXorf58, CXorf58

Immunogen:

Recombinant Human Putative uncharacterized protein CXorf58 protein (101-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A