

GPA33 Antibody, Biotin conjugated



PACO40205

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in cell-cell recognition and signaling.

Gene ID:

GPA33

Uniprot

Q99795

Synonyms:

Cell surface A33 antigen (Glycoprotein A33), GPA33

Immunogen:

Recombinant Human Cell surface A33 antigen protein (22-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A