

# MYCBP Antibody, HRP conjugated



PACO40219

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May control the transcriptional activity of MYC. Stimulates the activation of E box-dependent transcription by MYC.

**Gene ID:**

MYCBP

**Uniprot**

Q99417

**Synonyms:**

c-Myc-binding protein (Associate of Myc 1) (AMY-1), MYCBP, AMY1

**Immunogen:**

Recombinant Human c-Myc-binding protein (2-103AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A