

C12orf57 Antibody



PACO40246

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

In brain, may be required for corpus callusum development.

Gene ID:

C12orf57

Uniprot

Q99622

Synonyms:

Protein C10, C12orf57, C10

Immunogen:

Recombinant Human Protein C10 protein (2-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A