

SPIN2A Antibody, Biotin conjugated



PACO40269

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in the regulation of cell cycle progression.

Gene ID:

SPIN2A

Uniprot

Q99865

Synonyms:

Spindlin-2A (Protein DXF34) (Spindlin-like protein 2A) (SPIN-2) (SPIN-2A), SPIN2A, DXF34 SPIN2

Immunogen:

Recombinant Human Spindlin-2A protein (1-258AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A