

PACO40270

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in embryo hatching and implantation.

Gene ID:

Prss29

Uniprot

Q99MS4

Synonyms:

Serine protease 29 (EC 3.4.21. -) (Implantation serine proteinase 2) (ISP-2) (Strypsin-2) (Strypsin-related protein) (Tryptase-like proteinase), Prss29, lsp2

Immunogen:

Recombinant Mouse Serine protease 29 protein (18-179AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A