

# Prss29 Antibody, HRP conjugated



PACO40271

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in embryo hatching and implantation.

**Gene ID:**

Prss29

**Uniprot**

Q99MS4

**Synonyms:**

Serine protease 29 (EC 3.4.21. -) (Implantation serine proteinase 2) (ISP-2) (Strypsin-2) (Strypsin-related protein) (Tryptase-like proteinase), Prss29, lsp2

**Immunogen:**

Recombinant Mouse Serine protease 29 protein (18-179AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A