

PACO40318

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the reversible conversion of 3-phosphohydroxypyruvate to phosphoserine and of 3-hydroxy-2-oxo-4-phosphonooxybutanoate to phosphohydroxythreonine.

**Gene ID:**

Psat1

**Uniprot**

Q99K85

**Synonyms:**

Phosphoserine aminotransferase (PSAT) (EC 2.6.1.52) (Endometrial progesterone-induced protein) (EPIP) (Phosphohydroxythreonine aminotransferase), Pstat1, Psa Pstat

**Immunogen:**

Recombinant Mouse Phosphoserine aminotransferase protein (1-370AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A