

PELO Antibody, FITC conjugated



PACO40340

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for normal chromosome segregation during cell division and genomic stability. May function in recognizing stalled ribosomes and triggering endonucleolytic cleavage of the mRNA, a mechanism to release non-functional ribosomes and degrade damaged mRNAs. May have ribonuclease activity (Potential).

Gene ID:

PELO

Uniprot

Q9BRX2

Synonyms:

Protein pelota homolog (EC 3.1. -. -), PELO

Immunogen:

Recombinant Human Protein pelota homolog protein (1-385AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A