

ATG10 Antibody, Biotin conjugated



PACO40397

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

E2-like enzyme involved in autophagy. Acts as an E2-like enzyme that catalyzes the conjugation of ATG12 to ATG5. ATG12 conjugation to ATG5 is required for autophagy. Likely serves as an ATG5-recognition molecule. Not involved in ATG12 conjugation to ATG3. Plays a role in adenovirus-mediated cell lysis.

Gene ID:

ATG10

Uniprot

Q9H0Y0

Synonyms:

Ubiquitin-like-conjugating enzyme ATG10 (EC 2.3.2. -) (Autophagy-related protein 10) (APG10-like), ATG10, APG10L

Immunogen:

Recombinant Human Ubiquitin-like-conjugating enzyme ATG10 protein (1-50AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A