

C11orf68 Antibody, FITC conjugated



PACO40400

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

poly(A) RNA binding.

Gene ID:

C11orf68

Uniprot

Q9H3H3

Synonyms:

UPF0696 protein C11orf68 (Basophilic leukemia-expressed protein Bles03) (Protein p5326), C11orf68, BLES03

Immunogen:

Recombinant Human UPF0696 protein C11orf68 protein (1-251AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A