

# ZSCAN16 Antibody



PACO40410

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZSCAN16

**Uniprot**

Q9H4T2

**Synonyms:**

Zinc finger and SCAN domain-containing protein 16 (Zinc finger protein 392) (Zinc finger protein 435), ZSCAN16, ZNF392 ZNF435

**Immunogen:**

Recombinant Human Zinc finger and SCAN domain-containing protein 16 protein (1-348AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A