

# MRPL20 Antibody



PACO40510

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:****Gene ID:**

MRPL20

**Uniprot**

Q9BYC9

**Synonyms:**

39S ribosomal protein L20, mitochondrial (L20mt) (MRP-L20) (Mitochondrial large ribosomal subunit protein bL20m), MRPL20

**Immunogen:**


Recombinant Human 39S ribosomal protein L20, mitochondrial protein (46-149AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human skin tissue using PACO40510 at dilution of 1:100.