

MRPL20 Antibody, FITC conjugated



PACO40512

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

MRPL20

Uniprot

Q9BYC9

Synonyms:

39S ribosomal protein L20, mitochondrial (L20mt) (MRP-L20) (Mitochondrial large ribosomal subunit protein bL20m), MRPL20

Immunogen:

Recombinant Human 39S ribosomal protein L20, mitochondrial protein (46-149AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A