

# PPIL3 Antibody, HRP conjugated



PACO40555

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. May be involved in pre-mRNA splicing.

**Gene ID:**

PPIL3

**Uniprot**

Q9H2H8

**Synonyms:**

Peptidyl-prolyl cis-trans isomerase-like 3 (PPIase) (EC 5.2.1.8) (Cyclophilin J) (CyPJ) (Cyclophilin-like protein PPIL3) (Rotamase PPIL3), PPIL3

**Immunogen:**

Recombinant Human Peptidyl-prolyl cis-trans isomerase-like 3 protein (2-161AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A