

IFT22 Antibody



PACO40594

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Small GTPase-like component of the intraflagellar transport (IFT) complex B.

Gene ID:

IFT22

Uniprot

Q9H7X7

Synonyms:

Intraflagellar transport protein 22 homolog (Rab-like protein 5), IFT22, RABL5


Immunogen:

Recombinant Human Intraflagellar transport protein 22 homolog protein (1-185AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO40594 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human heart tissue using PACO40594 at dilution of 1:100.