

# ZFP64 Antibody



PACO40634

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

May be involved in transcriptional regulation.

**Gene ID:**

ZFP64

**Uniprot**

Q9NPA5

**Synonyms:**

Zinc finger protein 64 homolog, isoforms 1 and 2 (Zfp-64) (Zinc finger protein 338), ZFP64, ZNF338

**Immunogen:**


Recombinant Human Zinc finger protein 64 homolog, isoforms 1 and 2 protein (70-342AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Mouse liver tissue, Mouse kidney tissue. All lanes: ZFP64 antibody at 4 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 75, 69, 73, 49, 47 kDa. Observed band size: 75 kDa.