

PACO40674

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Regulates HCFC1 activity by modulating its subcellular localization. Overexpression of HCFC1R1 leads to accumulation of HCFC1 in the cytoplasm. HCFC1R1-mediated export may provide the pool of cytoplasmic HCFC1 required for import of virion-derived VP16 into the nucleus.

**Gene ID:**

HCFC1R1

**Uniprot**

Q9NWW0

**Synonyms:**

Host cell factor C1 regulator 1 (HCF-1 beta-propeller-interacting protein), HCFC1R1, HPIP

**Immunogen:**


Recombinant Human Host cell factor C1 regulator 1 protein (1-119AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung cancer using PACO40674 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO40674 at dilution of 1:100.