

CCDC94 Antibody



PACO40730

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

CCDC94

Uniprot

Q9BW85

Synonyms:

Coiled-coil domain-containing protein 94, CCDC94


Immunogen:

Recombinant Human Coiled-coil domain-containing protein 94 protein (1-279AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO40730 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO40730 at dilution of 1:100.