

Gkn2 Antibody



PACO40766

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

basal part of cell.

Gene ID:

Gkn2

Uniprot

Q9CQS6

Synonyms:

Gastrokine-2 (Blottin), Gkn2, Blot

Immunogen:

Recombinant Mouse Gastrokine-2 protein (21-184AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A