

# LHPP Antibody, FITC conjugated



PACO40776

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Phosphatase that hydrolyzes imidodiphosphate, 3-phosphohistidine and 6-phospholysine. Has broad substrate specificity and can also hydrolyze inorganic diphosphate, but with lower efficiency.

**Gene ID:**

LHPP

**Uniprot**

Q9H008

**Synonyms:**

Phospholysine phosphohistidine inorganic pyrophosphate phosphatase (hLHPP) (EC 3.1.3. -) (EC 3.6.1.1), LHPP

**Immunogen:**

Recombinant Human Phospholysine phosphohistidine inorganic pyrophosphate phosphatase protein (1-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A