

# C12orf49 Antibody, Biotin conjugated



PACO40817

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

C12orf49

**Uniprot**

Q9H741

**Synonyms:**

UPF0454 protein C12orf49, C12orf49

**Immunogen:**

Recombinant Human UPF0454 protein C12orf49 protein (34-205AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A