

ZMAT4 Antibody



PACO40818

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

p53 binding, RNA binding, intrinsic apoptotic signaling pathway by p53 class mediator.

Gene ID:

ZMAT4

Uniprot

Q9H898

Synonyms:

Zinc finger matrin-type protein 4, ZMAT4


Immunogen:

Recombinant Human Zinc finger matrin-type protein 4 protein (1-153AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO40818 at dilution of 1:100.