

NRN1 Antibody



PACO40846

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Promotes neurite outgrowth and especially branching of neuritic processes in primary hippocampal and cortical cells.

Gene ID:

NRN1

Uniprot

Q9NPD7

Synonyms:

Neuritin, NRN1, NRN


Immunogen:

Recombinant Human Neuritin protein (28-116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO40846 at dilution of 1:100.