

# TRIM36 Antibody



PACO40862

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

E3 ubiquitin-protein ligase which mediates ubiquitination and subsequent proteasomal degradation of target proteins. Involved in chromosome segregation and cell cycle regulation. May play a role in the acrosome reaction and fertilization.

**Gene ID:**

TRIM36

**Uniprot**

Q9NQ86

**Synonyms:**

E3 ubiquitin-protein ligase TRIM36 (EC 2.3.2.27) (RING finger protein 98) (RING-type E3 ubiquitin transferase TRIM36) (Tripartite motif-containing protein 36) (Zinc-binding protein Rbcc728), TRIM36, RBCC728 RNF98

**Immunogen:**


Recombinant Human E3 ubiquitin-protein ligase TRIM36 protein (70-280AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO40862 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO40862 at dilution of 1:100.