

# TAS2R14 Antibody




PACO40914

---

<b>Size:</b>	<b>Gene ID:</b>
50ug	TAS2R14
<b>Reactivity:</b>	<b>Uniprot</b>
Human	Q9NYV8
<b>Source:</b>	<b>Synonyms:</b>
Rabbit	Taste receptor type 2 member 14 (T2R14) (Taste receptor family B member 1) (TRB1), TAS2R14
<b>Isotype:</b>	<b>Immunogen:</b>
IgG	Recombinant Human Taste receptor type 2 member 14 protein (283-317AA).
<b>Applications:</b>	<b>Immunogen Species:</b>
ELISA, IHC	Human
<b>Recommended dilutions:</b>	<b>Storage: Buffer:</b>
IHC:1:20-1:200	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Form:</b>	<b>Storage:</b>
Liquid	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Conjugate:</b>	<b>Purification Method:</b>
Non-conjugated	>95%, Protein G purified
<b>Clonality:</b>	
Polyclonal	
<b>Research Area:</b>	
Neuroscience	

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO40914 at dilution of 1:100.