

# PDP2 Antibody, HRP conjugated



PACO40939

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the dephosphorylation and concomitant reactivation of the alpha subunit of the E1 component of the pyruvate dehydrogenase complex.

**Gene ID:**

PDP2

**Uniprot**

Q9P2J9

**Synonyms:**

[Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 2, mitochondrial (PDP 2) (EC 3.1.3.43) (Pyruvate dehydrogenase phosphatase catalytic subunit 2) (PDPC 2), PDP2, KIAA1348

**Immunogen:**

Recombinant Human [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 2, mitochondrial protein (67-529AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A