

PDP2 Antibody, Biotin conjugated



PACO40941

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the dephosphorylation and concomitant reactivation of the alpha subunit of the E1 component of the pyruvate dehydrogenase complex.

Gene ID:

PDP2

Uniprot

Q9P2J9

Synonyms:

[Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 2, mitochondrial (PDP 2) (EC 3.1.3.43) (Pyruvate dehydrogenase phosphatase catalytic subunit 2) (PDPC 2), PDP2, KIAA1348

Immunogen:

Recombinant Human [Pyruvate dehydrogenase [acetyl-transferring]]-phosphatase 2, mitochondrial protein (67-529AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A