

XPR1 Antibody, Biotin conjugated



PACO40945

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in phosphate homeostasis. Mediates phosphate export from the cell. May function in G-protein coupled signal transduction.

Gene ID:

XPR1

Uniprot

Q9UBH6

Synonyms:

Xenotropic and polytropic retrovirus receptor 1 (Protein SYG1 homolog) (Xenotropic and polytropic murine leukemia virus receptor X3) (X-receptor), XPR1, SYG1 XR

Immunogen:

Recombinant Human Xenotropic and polytropic retrovirus receptor 1 protein (529-696AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A