

LRRC61 Antibody



PACO40974

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

LRRC61

Uniprot

Q9BV99

Synonyms:

Leucine-rich repeat-containing protein 61, LRRC61


Immunogen:

Recombinant Human Leucine-rich repeat-containing protein 61 protein (1-259AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO40974 at dilution of 1:100.