

# TEKT3 Antibody



PACO40990

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules. Required for progressive sperm mobility.

**Gene ID:**

TEKT3

**Uniprot**

Q9BXF9

**Synonyms:**

Tektin-3, TEKT3

**Immunogen:**


Recombinant Human Tektin-3 protein (201-400AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver tissue using PACO40990 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO40990 at dilution of 1:100.