

IDI2 Antibody, HRP conjugated



PACO40999

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the 1,3-allylic rearrangement of the homoallylic substrate isopentenyl (IPP) to its highly electrophilic allylic isomer, dimethylallyl diphosphate (DMAPP).

Gene ID:

IDI2

Uniprot

Q9BXS1

Synonyms:

Isopentenyl-diphosphate delta-isomerase 2 (EC 5.3.3.2) (Isopentenyl pyrophosphate isomerase 2) (IPP isomerase 2) (IPPI2), IDI2

Immunogen:

Recombinant Human Isopentenyl-diphosphate delta-isomerase 2 protein (1-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A