

# IDI2 Antibody, Biotin conjugated



PACO41001

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the 1,3-allylic rearrangement of the homoallylic substrate isopentenyl (IPP) to its highly electrophilic allylic isomer, dimethylallyl diphosphate (DMAPP).

**Gene ID:**

IDI2

**Uniprot**

Q9BXS1

**Synonyms:**

Isopentenyl-diphosphate delta-isomerase 2 (EC 5.3.3.2) (Isopentenyl pyrophosphate isomerase 2) (IPP isomerase 2) (IPPI2), IDI2

**Immunogen:**

Recombinant Human Isopentenyl-diphosphate delta-isomerase 2 protein (1-227AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A